

Laboratory Liquid Scintillation Counting System

Model – SSS-22-LAB

Gamma Background Radiation Rejection Features

- Energy analyzer window rejects pulses with energies outside the window setting
- Optional Lead Shielding around detector

Counting Assembly Features

- Excellent repeatability
- Fully light tight system
- Fail safe interlock to protect PM tubes
- High transmission optical coupling to PM tube

PM Tube and Pre-Amp Noise Elimination Features

- Coincidence Verification with High quality Dual PM tubes and preamps.
- Fully adjustable energy analyzer window rejects low energy pulses.
- **SSS-22-LAB:** 3 Vial capacity for faster thru-put and easy comparison of sample to calibration standard or to the background – 2 PM Tubes
- 3 Vial capacity provides less wait time for dark adaption of phosphor.

Data Analysis and Presentation

- Scintillation counts detected by PM tubes are processed by a fully adjustable single channel analyzer which is centered on the energy peak of the isotope being measured. This deletes both higher energy pulses from background radiation and lower energy counts from the PM tube or circuit noise. Pulses feed to a digital scaler and optional digital printer and/or user's PC.
- This process allows long count times for measurement of very minute samples.
- R232 Interface with easy conversion to USB interface to computers or data stations.

SPECIFICATIONS:

H-3 Gross Efficiency:	>50% Gross, Coincidence efficiency >10% Verified.
Count Times:	1 sec. thru 99 sec. (1 sec. increments), and 1 min. thru 99 min. (in 1 min. increments).
Voltage:	0-1500 Volts - fully user settable.
Readout:	Digital - 6 digit LCD, (LED optional).
Outputs:	Standard: Serial pulse output and R-232.
Power:	Rechargeable Lithium batteries and AC adapter
Sample Size:	Accepts standard Liquid Scintillation vials 20-25 ml.
Case:	Weather-tight NEMA rated tough, durable, lightweight. A "tongue-in-groove" continuous "O"-ring system makes the case water, moisture, and dust tight
SSS-22-LAB-SS Case:	Stainless Steel Weather-tight, easy to clean
OPTIONAL:	SSS-22-LAB: Accepts 6-7 ml vials with Optional adapter.
Scintillation Fluors:	Accepts most scintillation fluors. Perkin-Elmer Ultima Gold-LLT is recommended for H-3 counting



TECHNICAL ASSOCIATES OVERHOFF TECHNOLOGY

7051 ETON AVENUE, CANOGA PARK, CALIFORNIA 91303
P: 818-883-7043 | C: 818-836-0381 | F: 818-883-6103
SALES@TECH-ASSOCIATES.COM | TECH-ASSOCIATES.COM

DIVISIONS OF



US NUCLEAR CORP

TICKER UCLE